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**Kaynaklar**

- Brisson-Noël A, Delrieu P, Samain D, Courvalin P. Inactivation of lincosaminide antibiotics in *Staphylococcus*. Identification of lincosaminide O-nucleotidyltransferases and comparison of the corresponding resistance genes. *J Biol Chem* 1988; 263(31):15880-7
- Chambers HF, Hackbarth CJ, Drake TA, Rusnak MG, Sande MA. Endocarditis due to methicillin-resistant *Staphylococcus aureus* in rabbits: expression of resistance to beta-lactam antibiotics in vivo and in vitro. *J Infect Dis* 1984; 149(6):894-903.
- Hombach M, Weissert C, Senn MM, Zbinden R. Comparison of phenotypic methods for the detection of penicillinase in *Staphylococcus aureus* and proposal of a practical diagnostic approach. *J Antimicrob Chemother* 2017; 72(4):1089-1093. DOI: 10.1093/jac/dkw521.
- Kaatz GW, Seo SM. Mechanisms of fluoroquinolone resistance in genetically related strains of *Staphylococcus aureus*. *Antimicrob Agents Chemother* 1997; 41(12):2733-7.
- LaPlante KL, Leonard SN, Andes DR, Craig WA, Rybak MJ. Activities of clindamycin, daptomycin, doxycycline, linezolid, trimethoprim-sulfamethoxazole, and vancomycin against community-associated methicillin-resistant *Staphylococcus aureus* with inducible clindamycin resistance in murine thigh infection and in vitro pharmacodynamic models. *Antimicrob Agents Chemother* 2008; 52(6):2156-62. DOI: 10.1128/AAC.01046-07.
- Mendes RE, Farrell DJ, Flamm RK, Sader HS, Jones RN. Analysis of Vancomycin Susceptibility Testing Results for Presumptive Categorization of Telavancin. *J Clin Microbiol* 2015; 53(8):2727-30. DOI: 10.1128/JCM.00611-15.
- Papanicolas LE, Bell JM, Bastian I. Performance of phenotypic tests for detection of penicillinase in *Staphylococcus aureus* isolates from Australia. *J Clin Microbiol* 2014; 52(4):1136-8. DOI: 10.1128/JCM.03068-13.
- Peñuelas M, Candel FJ, Lejarraga C, López-González L, Viñuela-Prieto JM, López de Mendoza D. Activity of linezolid and tedizolid against clinical isolates of methicillin-resistant and methicillin and linezolid resistant *Staphylococcus aureus*: an in vitro comparison. *Rev Esp Quimioter* 2016; 29(5):255-8
- Sierra JM, Cabeza JG, Ruiz Chaler M, Montero T, Hernandez J, Mensa J, Llagostera M, Vila J. The selection of resistance to and the mutagenicity of different fluoroquinolones in *Staphylococcus aureus* and *Streptococcus pneumoniae*. *Clin Microbiol Infect* 2005; 11(9):750-8.
- Skov R, Larsen AR, Kearns A, Holmes M, Teale C, Edwards G, Hill R. Phenotypic detection of mecC-MRSA: cefoxitin is more reliable than oxacillin. *J Antimicrob Chemother* 2014 Jan;69(1):133-5. doi: 10.1093/jac/dkt341. Epub 2013 Sep 12.
- Trzcinski K, Cooper BS, Hryniewicz W, Dowson CG. Expression of resistance to tetracyclines in strains of methicillin-resistant *Staphylococcus aureus*. *J Antimicrob Chemother* 2000; 45(6):763-70.
-